

Chromai

Innovation for Better Life

Website: www.chromai.com

Chromai · Beijing

Address: 5th Floor, Building 10, Liandong U Valley · Lanbei Technology Park, No. 3 Xixing Road,
Shunyi District, Beijing, China

Tel: (86-10)80429960

Chromai · Suzhou

Address: 5th Floor, Building 5, No. 1268 Jiaotong Road, Wujiang District, Suzhou, Jiangsu, China

Tel: (86-512)63226066

Email: contact@chromai.com



WeChat
Official Account



Linkedin
Official Account

Product information in this document may change without notice. For the latest updates on Chromai's products, applications, and services, please visit the official website or contact us.

Chromai®

— 科诺美 —

UHPLC

Ultra-High Performance
Liquid Chromatography Series

Innovation for Better Life

30%

R&D Technical Staff

10+

Product Series

10+

Application Fields

30+

Independent Intellectual
Properties

100+

Global Coverage
Areas

Chromai , as a leader in China's high-end scientific instruments and medical diagnostic equipment, is committed to making analytical testing more efficient, precise, and user-friendly. With a focus on independent R&D and continuous innovation, the company possesses advanced technologies and numerous core patents. It specializes in integrating scientific instruments with automation and intelligence, offering products and solutions widely used by global clients across multiple fields. Chromai is a national high-tech enterprise and a member of organizations such as the China Instrument and Control Society, China Association for Instrumental Analysis, and China Medical Device Industry Association. The company strives to elevate China's manufacturing reputation worldwide, becoming a Chinese tech company that promotes global progress in health and safety while improving efficiency, reducing analytical costs, and creating value for customers.



Change the global industry market pattern, let China create a global reputation

Efficient and reliable solutions
Empowering laboratory analytical performance

Every data reliable and trust worthy

Advantages
and Value

Scene
adaptation

Adaptation through a full-scenario LC matrix — from analytical ultra-high-pressure systems to preparative purification platforms, supplemented with specialized modules to precisely address complex sample detection challenges, meeting demands from scientific research to production.

Improve
efficiency

Ultra-high-pressure separation cuts analysis time by 60%. High-sensitivity detection breaks quantitation limits.

Service
Value-Added

3Q certification and 24/7 support. Optimizes asset allocation, streamlines upgrades, and ensures a balance between research efficiency and cost control—maximizing ROI.

Diverse Scenarios
and Customized Solutions

- **Environmental Monitoring**
Pollutant Monitoring
- **Agriculture & Livestock**
Pesticides | Veterinary Drugs | Feed
- **Research & Education**
University Labs | Research Institutes
- **Food Safety**
Nutrients | Harmful Substances |
Illegal Additives
- **Biopharma & Pharmaceuticals**
Biologics | Chemical Drugs | Synthetic Biology |
Traditional Medicine | Excipients
- **Chemical Materials**
Cosmetics | Polymers | Fuels
- **Regulatory Testing**
FDA | Customs | Agricultural Inspection |
Quality Inspection | Commercial Inspection
- **Forensics & Judiciary**
Judicial Identification | Forensic Toxicology
- **Public Health**
Disease Control

Chromai
— 科诺美 —



Leaps
UHPLC System



Leaps Pro
UHPLC System Pro



Frontier
UHPLC System



Leaps 2D
2D UHPLC System



Leaps Pro Bio
Bio-Inert LC System



Leaps Prep
Preparative LC System

Comprehensive and powerful UHPLC product lineup, from rapid analysis to complex samples, bio-inert to preparative purification.
Full coverage of column particle size selection: 1.7 μ m for UHPLC, 3~5 μ m for routine analysis, 7~20 μ m for preparative purification.

Multi-Domain

Multi-Module

Reliable and
Stable

Cost-Effective

Analytical

Preparative

Rapid Analysis with Reliable Performance

10500PSI (720 Bar)

Maximum Pressure

2x54 / 2x96 Well Plates

Sample Tray

Injection synchronization

Improves retention time accuracy

Leaps



Dual ternary pumps

Multifunctional modules

Switch between standard LC and 2D LC

Supports custom solutions

Online enrichment and impurity removal

Minimizes manual sample preparation

Leaps 2D

● Ultra-stable pulse control for precise operation

Designed for high stability—maintains accurate performance even with low-ratio mixing; pressure pulsation is controlled within <10 PSI to ensure reliable operation

● Fully automated sample processing

Integrates pre-column derivatization, automatic dilution, and mixing—enabling full-process automation

● Flexible and adjustable multi-mode configuration

Equipped with a large-capacity column oven that holds up to 4x300mm columns; temperature control range of 4°C to 85°C; supports 10-/6-port valves for multi-center cut functionalities

● Fully intelligent adaptive operation

Supports cross-sequence operations—automatically analyzes multiple projects in sequence without manual

Top LC Sytem in China

10500 PSI (720 Bar)

Maximum Pressure

4×54 / 4×96 Well Plates

Sample Tray

Hand-Tight Zero Dead
Volume Tubing

Minimize Extra-Column Volume

Leaps Pro



22040 PSI (1520 Bar)

Maximum pressure

Total reflection flow cell

Low Dispersion

LC-MS Combinations

SDK System Integration

Frontier

High-efficiency analytical experience

Smart delay volume reduction tech, combined with zero dead volume tubing and semi-micro flow cells to Minimize Extra-Column Volume

Smaller volume and Non-Destructive Sampling

Integrated Loop Autosampler, minimum injection volume to 0.1 μ L, supports multi-solvent washing, reduces sample carryover to 0.0025%

Precise temperature control system

Semiconductor temp. control with a range of 4°C – 85°C, built-in high-sensitivity leak sensor, optional preheating module available

Compatible with multiple detector types

Support Various Detectors such as UVD, DAD, FLD, ELSD*, and RID*

Leaps Pro Bio

Bio-inert LC System


Designed for efficient analysis of biological samples

- Inert flow path for biological samples
- Biocompatible materials
- Compatible with biological samples and extreme pH conditions
- Reduced metal adsorption




Fully inert and corrosion-resistant flow path: Accurate analysis with hassle-free operation

- Flow path with corrosion-resistant materials, enabling long-term stability in high-salt or extreme pH
- Reducing nonspecific adsorption and sample loss, without complex mobile phases or frequent passivation
- Fully inert configuration including pump, flow path, injection valve, syringe, and flow cell
- Plunger wash peristaltic pump to effectively reduce salt precipitation
- Hand-tight zero-dead-volume fittings prevent diffusion and peak broadening



Vairous application scenarios and analytical modes

- Supports RP, IEX, SEC, and HIC modes
- mAb/BsAb purity analysis, mRNA impurity detection, and quality control of CGT



High-precision automated sample handling

- Customized with high-precision robotic arm modules integrated with recognition and path-planning algorithms, enabling automated sample tray replacement and precise slot positioning for full end-to-end unmanned operation



Leaps Prep

Preparative LC System
Simple and flexible preparative purification

- Fully automated workflow
- Flexible expandability
- Stable separation



High-pressure preparative applications


- Flow Rate: 0–150 mL/min
- Pressure: 0–400 Bar
- Supports up to 4 solvents for versatile application scenarios

Enables fully automated fraction collection

- Automated sample injection and fraction collection
- Various collection modes: Time, Threshold, Slope, and Mixed Mode
- Various rack types, including 30 mm, 16 mm, 13 mm Tubes, and 96-well plates

Large-volume injection and column switching

- Multiple injection modes, with a maximum injection volume of 5 mL
- Optional refrigerated sample compartment available
- Loading Ranges from Milligrams to Grams
- Flexible sample tray options: 2/5/7 mL sample vials
- Upgradeable to 2D Preparative LC System




High-pressure preparative pump




● Diode Array Detector (DAD)



- 1024 pixels, wavelength range of 190—800nm; reference wavelength eliminates background interference and drift correction function
- Slit widths of 1.2nm and 6nm for high-resolution spectra and more accurate monitoring data
- Simultaneous monitoring of 12-channel wavelengths
- 3D spectral monitoring and peak purity calculation function
- Supports exporting chromatograms at any wavelength within the spectral scanning range



- Ultra-high sensitivity, low dispersion total reflection flow cell, Designed to ultra-high efficiency analysis
- Maximum path length of 60mm for 6x signal sensitivity



Comparison of 10 mm and 60 mm total reflection flow cells


● UV-Vis Detector




- Wavelength range of 190—900nm
- 100Hz data rate enables accurate recovery of narrow peaks
- 2.5AU linear range
- Wavelength-time programming, wavelength switching, and stop-flow scanning
- Light source is easy to replace, various flow cell optional

● FLD

- New optical design
- Delivers excellent signal-to-noise ratio
- Wavelength range: 200—750nm/900nm
- Max.data rate: 100Hz



Aflatoxin standard solution



Polycyclic aromatic hydrocarbons (PAHs) standard solution



- **Food Safety:** Detecting harmful substances in food
- **Chemical Analysis:** Organic synthesis and pharmaceutical R&D
- **Life Sciences:** Glycosylation analysis and other applications
- **Environmental Monitoring:** Pollutant analysis in environmental samples

● Evaporative Light Scattering Detector (ELSD)

- Universal detector, suitable for gradient elution
- High sensitivity, low sensitivity to temperature changes, stable baseline
- Online filtration with adjustable pressure valve design, clean gas source and stable airflow
- Precise temperature control system
- Front-positioned nebulization chamber for easy cleaning
- Optional GPC module




Figure 1: Astragaloside IV standard




Figure 2: Ejiao Tonic Paste sample

Application Scope

Traditional Chinese Medicine

Glycosides / Alkaloids
Terpenoids / Flavonoids
/ Steroidal compounds

Chemical Pharmaceuticals

Aminoglycosides

Carbohydrates

Monosaccharides
/ Oligosaccharides
/ Sugar alcohols


Polymeric Surfactants

Polysorbate
/ Polyvinyl alcohol

Lipids


Phospholipids
/ Triglycerides
/ Fatty acids

● Refractive Index Detector (RID)



- High sensitivity for sugar analysis; detection limit up to 10^{-8} g/mL
- Temperature control and automatic flush-balance to quickly stabilize the baseline
- Optional GPC module

● Lotus Series Columns



- Packing material from globally reliable suppliers
- Industry-leading GMP standards, Meticulous column packing process
- Special functional group bonding density, Balances for hydrophobic and hydrophilic compounds
- Various application requirements, Particle sizes from 1.7μm to 20μm
- Leader in Bioinert LC Systems in China, Complete solutions for LC consumables
- Complex sample separation, High surface area and high carbon loading



EyouLab Intelligent Data System

Enhances laboratory efficiency with greater accuracy and speed
Effectively reduces operational costs and risks



EyouLab CDS

EyouLab fully complies with FDA 21 CFR Part 11, GMP, GLP, and other regulatory requirements. It offers multiple deployment options, including standalone, networked, and medical versions, with multilingual support covering everything from sample preparation to precision medicine. Built on dozens of proprietary technologies, it combines a “smart guided interface + cloud-based data management” to simplify complex experimental operations, ensure data integrity and accuracy, and significantly boost lab capabilities and efficiency—accelerating the transition to intelligent laboratories.



More user-friendly and intelligently adaptable for analytical workflows

- › Integrated normalization window with a minimalist interface simplifies complexity for ease of use
- › Graphical monitoring interface provides real-time status of instrument operation
- › Easy-to-use report editor with readily accessible templates
- › Robust data processing with over 34 integration functions and calculation methods
- › Advanced automatic peak width and threshold algorithms for optimal integration




Standalone Operation

Remote Access

Archiving Management

Audit Trail

Risk Management



More Specific and Customized Solutions

Precise adaptation and performance upgrade

Tailored solutions based on sample characteristics

Reduce operation and maintenance costs, meet compliance standards



Amino Acid Analyzer

Pre-column derivatization via auto-sampler and post-column derivatization with ninhydrin; the fully automated amino acid analyzer offers efficient automation for primary and secondary amino acid derivatization in pharmaceutical R&D and food testing, enhancing analysis efficiency and accuracy



Online Pre-treatment

Streamlines sample pre-treatment process, doubling detection efficiency. The dual ternary pump configuration saves space and ensures easy operation. Highly integrated design significantly optimizes workflow while reducing both cost and time




GPC/SEC

Utilizes high-accuracy, high-performance SEC/GPC systems for precise analysis of polymers and macromolecules, meeting the full-range needs from molecular structure research to quality control



2D Preparative Chromatography

Combines powerful separation performance with a simplified operational experience. With dual-mode compatibility, the system easily adapts to various application scenarios—from method development to routine analysis—unlocking broader possibilities in sample analysis



Scientific Research and Education

› Tandem LC

Loop and tandem LC utilize parallel column regeneration technology to save time on ineffective column cleaning and re-equilibration, significantly improving sample throughput

› Parallel LC

Parallel LC enables cross-gradient operation of two completely different LC methods simultaneously, increasing productivity without expanding workspace. It allows rapid acquisition of sample data while minimizing the risk of degradation of unstable analytes

Environmental Monitoring

› Detection of Polycyclic Aromatic Hydrocarbons in Water

Online SPE column are used, and enriched compounds are transferred to the analytical flow path via switching valves for quantitative analysis

Chemical Materials

› Determination of Compound Content in Cosmetics

Simultaneous determination of α -arbutin, β -arbutin, deoxyarbutin, hydroquinone, and phenol in cosmetics, with detection and quantification limits meeting testing requirements.

Food Safety

› Fat-Soluble Vitamins in Food

Using 2D-LC with heart-cutting, the tedious steps of VD normal-phase purification are avoided. After sample preparation, a single injection allows simultaneous detection of vitamins A/ D/E

Agriculture

› Analysis of Aflatoxins in Feed

Photochemical derivatization is low-cost and avoids the use of highly concentrated additives that can damage instruments in chemical derivatization. Compared with pre-column derivatization, it is safer and less labor-intensive

Pharmaceuticals & Biopharmaceuticals


› Purity and Impurities Analysis of Nucleic Acid Drugs

Anion exchange chromatography is used to quickly separate sense and antisense strands and impurities in mRNA drugs. The results are reliable and meet the requirements of nucleic acid drugs for HPLC column systems

› Analysis of Small-Molecule Toxins and ADCs

IEC in mAb and ADC workflows ensures high sensitivity in analysis and detection, supporting precise quality control of ADC drugs


Chromai is deeply engaged in the validation, optimization, and innovative development of analytical methodologies across various fields, continuously delivering cutting-edge industry application solutions and customized services to help analytical laboratories remain competitive and achieve flexible expansion. Our business now comprehensively covers diverse fields including academic research, pharmaceuticals, biopharmaceuticals, agriculture and animal husbandry, chemical materials, and environmental monitoring. With full-scenario technical service capabilities, we provide precise support for the analytical needs of various industries.



Full Lifecycle Service Assurance: Supporting You Every Step of the Way

From Deployment to Advancement, Professional Operation & Maintenance Ensures Worry-Free Support

To ensure continuity in experimental workflows and data reliability, we build a full lifecycle service system for every delivered instrument. From the moment your order is confirmed, a dedicated engineer will be assigned to provide services including installation, commissioning, and standardized training. During the instrument's lifecycle, our engineers will provide regular preventive maintenance and fast, responsive support, ensuring smooth operation of your instrument and delivering efficient and worry-free assurance.



Make Analysis and Detection More Efficient, Accurate and Easy To Use